

## Comfort<sup>®</sup> Personal Cleansing Rinse-Free Shampoo Cap



### Packaging Information

#### Contents:

- One premoistened shampoo cap.

#### Instructions:

Tear open label. If desired, warm entire package in the microwave (1000W) **no longer than 30 seconds**. Caution: Test temperature of cap prior to use. Wait if excessively hot. Place cap on head and massage until hair is fully saturated (1-3 minutes). Remove, towel dry, comb hair.

### Specifications

#### Product Information

Product name:	Comfort Personal Cleansing Rinse-Free Shampoo Cap
Product code:	7409-X
Description	One premoistened cap.
Accreditation	ISO 13485:2003
Storage Information	Do not store above 25°C.
Expiration Date	None

#### Ingredients

Purified water (aqua)  
Isostearamidopropyl morpholine lactate  
Disodium cocoamphodictetate  
DMDM hydantoin  
Phenoxyethanol  
Methylparaben  
Butylparaben  
Ethylparaben  
Propylparaben  
Fragrance (parfum)  
Simethicone

#### Product Specifications

Diameter of cap	27.94 cm
Material composition of outer layer	Polyethylene
Material composition of inner layer	Polyester
Amount of solution in cap	0.3625 lbs/5.8 oz/170mL
Amount of cap in a case	40